

# dtintegrate dtprofit

SBC01-????-00

January 9, 1996

Prepared by: J.W. Pflugrath  
Molecular Structure Corporation  
3200 Research Forest Drive  
The Woodlands, Texas 77381  
(713) 363-1033  
(713) 364-3628 FAX  
jwp@msc.com

Prepared for: Contract No. 943072401  
Mary Westbrook  
Argonne National Laboratory  
9700 South Cass Avenue  
Argonne, Illinois 60439  
(708) 252-8914  
(708) 252-4021 FAX  
westbrook@anlel.el.anl.gov

Reviewed by:

Approved by:

Copyright © 1996, 1995, 1994 Molecular Structure Corporation  
RESTRICTED RIGHT NOTICE SHORT FORM (JUNE 1987)

Use, reproduction, or disclosure is subject to restrictions set forth in Contract No. W-31  
109-ENG-38 and Contract No. 943072401 with the University of Chicago, Operator of  
Argonne National Laboratory

## TABLE OF CONTENTS

<b>1 SCOPE AND PURPOSE</b>	<b>3</b>
<b>2 BACKGROUND</b>	<b>3</b>
<b>3 INITIAL INPUT IMAGE</b>	<b>5</b>
<b>4 RUNNING DTINTEGRATE</b>	<b>7</b>
4.1. Command line options	7
4.2. Examples	9
4.3. Results	9
4.4. Overcoming difficulties	19
<b>5 FILE FORMATS</b>	<b>20</b>

## 1 Scope and Purpose

This document describes how to use the `dtintegrate` and `dtprofit` modules of d\*TREK. Together `dtintegrate` and `dtprofit` integrate a scan or series of single crystal X-ray diffraction images such as those created by the `dtcollect` module of d\*TREK to produce a reflection list that contains the Miller index (hkl), estimated intensity and standard deviation for the Bragg reflections that appear in the scan. `Dtintegrate` uses the same procedures as `dtpredict` to predict reflection positions and as `tCrefine` to refine crystal, detector and source properties. The documentation for these modules should be consulted for more information.

`dtintegrate` reads the individual images of a scan and collects three-dimensional profiles from predicted reflection positions. Every few images (set by the `-batch` command line argument), the differences between observed and calculated reflection positions are minimized. The reflection profiles are written to a file named `reflnprofiles.profit`. `dtprofit` reads this file and performs 3D profile analysis with the method of Kabsch (1988) and writes a reflection list file called `dtprofit.ref`.

`dtintegrate` gets its initial information about the detector, source and crystal properties from the header of an image file. This file is usually the output of the `tCrefine` module. This information is needed to predict reflection positions. The default is to integrate all reflections predicted in the scan, so that additional input is not required for integration.

### 1.1 Definitions and Abbreviations

The definitions described in the `dtcollect`, `tCrefine`, and `dtpredict` documents are appropriate for `dtintegrate`.

## 2 Background

In order to integrate the intensities of Bragg reflections in a single crystal diffraction experiment, their positions in images must be predicted accurately. This can be done provided information about the source, the crystal, the crystal goniometer, the detector, the detector goniometer, and the rotation is known. For the predictions to be accurate, the crystal, source and detector properties must be refined by minimizing the differences between observed and predicted reflection positions. In order to accurately predict reflections for images collected after the start of a scan, properties of the crystal, detector, and source should be refined periodically throughout integration. One can assume the detector and source properties do not change over the course of a scan data collection and therefore not refine them. However, experience shows that the initial values for the

detector and source may not be the most valid, so it is best to refine them just as one would refine crystal properties. Crystal properties can change due to radiation damage and incorrectly mounted crystals (i.e., the crystal can dry out or slip which changes the unit cell parameters and the crystal orientation).

Integration thus proceeds by looping through the images in a scan. For each image:

- A. If this image is the last one in a 'batch' of images, then refine the crystal, detector and source properties based on the observed positions of reflections.
- B. The reflections which appear in that image are predicted and the positions of the full reflections are re-predicted.
- C. For each predicted reflection:
  1. If this is the first time this reflection is predicted, then allocate a 3D array to hold a set of pixels that would completely encompass the reflection. The array should be large enough to include pixels for a local background determination and allow for small shifts of the predicted reflection position due to crystal slippage and inaccurate prediction.
  2. Next copy a window of pixels around the reflection from the image to the 3D array and check if the 3D array is completed.
  3. If the array is full, then pre-process the reflection. It is corrected for non-uniformity of response and the observed reflection centroid is determined if the peak is strong enough. Mark the reflection as full.
  4. Then go on to the next reflection.
- D. If this image is the last one in a 'batch' of images, then process the 'full' reflections. This means determine a local background, decide in a preliminary analysis which pixels are in background regions, which are in the peak region and which are in neither region (i.e. belong to the peak of another reflection). Process the array into 3D profiles. Write these profiles to disk for subsequent profile analysis.
- E. Then go on to the next image.

When all images have been read in and processed, the result is a reflection list with preliminary integrated intensities and a file containing reflection profiles. `dtprofit` is then run to do the actual profile fitting. In profile fitting, one makes the valid assumption that reflections have the same spot boundary and the same (except for a constant scale factor) intensity distribution within the boundary. In other words, reflections have the same expected profiles except each is scaled by a different constant. Thus if a reference profile is determined, then a least-squares fit between the observed profile and the reference profile can be made and the scale factor determined. If the reference profile is scaled such that its integrated intensity is 1, then the scale factor between observed and reference profile is equal to the estimated intensity of the reflection. For more information on profile fitting please consult the following references:

Diamond, R. (1969) *Acta Cryst.* **A25**, 43-55. Profile Analysis in Single Crystal Diffractometry.

Ford, G.C. (1974) *J. Appl. Cryst.* **7**, 555-564. Intensity Determination by Profile Fitting Applied to Precession Photographs.

Rossmann, M.G. (1979) *J. Appl. Cryst.* **12**, 225-238. Processing Oscillation Diffraction Data for Very Large Unit Cells with an Automatic Convolution Technique and Profile Fitting.

Rossmann, M.G., Leslie, A.G.W., Abdel-Meguid, S.S. and Tsukihara, T. (1979) *J. Appl. Cryst.* **12**, 570-581. Processing and Post-Refinement of Oscillation Camera Data.

Kabsch, W. (1988) *J. Appl. Cryst.* **21**, 916-924. Evaluation of Single-Crystal X-ray Diffraction Data from a Position-Sensitive Detector.

Observed intensities are corrected for the Lorentz and polarization factors calculated by the prediction algorithm. They are not otherwise scaled, however. A subsequent step is needed to merge, scale and average the reflections before they are used in any crystallographic calculations.

### 3 Initial input image

`dtintegrate` reads the header of an input image specified on the command line to get the initial values for most of the objects that it manipulates. These include the source, crystal, detector, crystal goniometer, and scan properties. In order to begin its work, `dtintegrate` does not read the binary pixel data in the input image which may be even be missing (i.e. the image is dimensioned 0 by 0). If the detector spatial distortion and non-uniformity of response information in a header requires additional images to be read in, then `dtintegrate` does so automatically. A typical image header is shown below. The keywords required by `dtintegrate` are the same as those required by `dtpredict`, so please refer to that documentation for their descriptions. The close relationship between `dtintegrate` and `dtpredict` is obvious once you realize that the routines of `dtpredict` are used by `dtintegrate` to predict reflections.

```
1 {
2  HEADER_BYTES= 2048;
3  TYPE=mad;
4  SIZE1=0;
5  SIZE2=0;
6  CRYSTAL_GONIO_VALUES=0.0 0.0 0.0;
7  COMMENT=Header edited by dtheadredit;
8  SCAN_TITLE=start end inc time nOsc nDark nDup nDlim nDC nDCup;
1  SCAN_ROTATION=0.0 12.0 0.2 4 0 1 0 100 1 0;
2  SCAN_ROTATION_AXIS_NAME=omega;
3  SCAN_ROTATION_VECTOR=1.0 0.0 0.0;
4  SCAN_TEMPLATE=lyso3.####.ss;
```

```
5  SCAN_SEQ_INFO=0 2 0;
6  ROTATION=0.0 0.2 0.2 4 0 1 0 100 1 0;
7  ROTATION_VECTOR=1.0 0.0 0.0;
8  ROTATION_AXIS_NAME=omega;
9  CRYSTAL_UNIT_CELL=88.00 88.00 104. 90.00 90.0 120.0;
10 CRYSTAL_MOSAICSPREAD=0.145;
11 CRYSTAL_DESCRIPTION=Test of dtpredict;
12 CRYSTAL_SPACEGROUP=175;
13 CRYSTAL_REFINE_FLAGS=0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0;
14 CRYSTAL_ORIENT_ANGLES=152.3 5.0 152.9;
15 CRYSTAL_ORIENT_VECTORS=1 0 0 0 1 0 0 0 1;
16 SOURCE_VECTORS=0.0 0.0 1.0 0 1 0 1 0 0;
17 SOURCE_POLARZ=0.5 1.0 0.0 0.0;
18 SOURCE_SPECTRAL_DISPERSION=0.0002 0.0002;
19 SOURCE_SIZE=0.0 0.0 0.0 0.0;
20 SOURCE_CROSSFIRE=0.0 0.0 0.0 0.0;
21 SOURCE_WAVELENGTH=1 1.54178;
22 SOURCE_REFINE_FLAGS=0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0;
23 DETECTOR_NUMBER=1;
24 DETECTOR_NAMES=D0_;
25 D0_NONUNF_TYPE=Dark_only;
26 D0_NONUNF_INFO=nonunf.img dark.img;
27 D0_SPATIAL_DISTORTION_TYPE=Interp_spatial;
28 D0_SPATIAL_DISTORTION_INFO=distor;
29 D0_DETECTOR_DIMENSIONS=512 512;
30 D0_DETECTOR_SIZE=50.0 50.0;
31 D0_DETECTOR_VECTORS=1 0 0 0 1 0;
32 D0_DETECTOR_DESCRIPTION=ANL-SBC gold detector single module;
33 D0_DETECTOR_REFINE_FLAGS=0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0;
34 D0_GONIO_NUM_VALUES=6;
35 D0_GONIO_NAMES=RotZ RotX RotY TransX TransY TransZ;
36 D0_GONIO_UNITS=deg deg deg mm mm mm;
37 D0_GONIO_VECTORS=0 0 1 1 0 0 0 1 0 1 0 0 0 1 0 0 0 -1;
38 D0_GONIO_VALUES=1.0 2.0 0.0 -.5 -.8 102.3;
39 D0_GONIO_DESCRIPTION=A simple detector goniostat;
40 CRYSTAL_GONIO_NUM_VALUES=3;
41 CRYSTAL_GONIO_NAMES=Omega Chi Phi;
42 CRYSTAL_GONIO_UNITS=deg deg deg;
43 CRYSTAL_GONIO_VECTORS=1 0 0 0 1 0 1 0 0;
44 CRYSTAL_GONIO_DESCRIPTION=A 3-circle eulerian crystal goniostat;
45 FILENAME=hex3.img;
46 BYTE_ORDER=big_endian;
47 Data_type=short int;
48 COMPRESSION=None;
49 }
```

Other command line options allow one set parameters specific to dtintegrate. These are described in the section 4.

## 4 Running dtintegrate

After dtintegrate has been installed and placed in your PATH, just enter dtintegrate along with command line options to run it. It parses the command line arguments, reads any required input images, spatial distortion files, non-uniformity files and writes a reflection profile file if reflections are integrated. Messages are written to stdout and stderr as required. The syntax for running dtintegrate is:

```
dtintegrate [-scan scanfile][ -reso rmin rmax |  
-batch num_images_per_batch |  
-seq start_seq_num end_seq_num |  
-pad num_pad_images  
-window fast_size slow_size  
-profsize profile_size1 profile_size2  
-nonunf file_with_nonunf_info_in_header  
]
```

dtprofit

### 4.1 Command line options

Although dtintegrate will run without any command line options if scan information is found in the header of a file called scan.img in the current working directory, most situations will benefit from the use of one or more command line options described below. dtprofit has no command line options.

#### Command line options: Description and default

-scan scanfile

The header of image scanfile (no default name) is read for the keywords described in section 3 to set up the experiment properties required for predicting reflections. The keywords SCAN\_ROTATION, SCAN\_ROTATION\_VECTOR are used to specify the rotation properties. If a -scan option is not specified, then dtintegrate uses '-scan scan.img' as the default.

-reso fReso1 fReso2

A resolution range in Ångstroms which overrides the default resolution range determined from the source wavelength and the detector(s) position. In order not miss any reflections on the image, do not use this option. If fReso1 is 0.0 or less, then the resolution range is also determined from the source wavelength and detector(s) position.

-seq Start\_seq\_num End\_seq\_num

This overrides the default sequence of images that will be integrated as determined from the SCAN\_ROTATION keyword. Start\_seq\_num is the

sequence number of the first image to be integrated while `End_seq_num` is the sequence number of the last image to be integrated. The default is to integrate every image specified by the `SCAN_ROTATION` keyword in the scan file. Use this command line option to choose a specific sequence of contiguous images in a scan.

`-batch num_images_per_batch |`

This overrides the number of images (default=5) that will be processed before a refinement is carried out and reflection 3D profiles generated. It wastes time to refine too often, unless the crystal properties are changing drastically. It is also unwise to have the batch size very large as the 3D reflection profiles are kept in core memory until after a refinement is done so that accurate centroids for them can be re-predicted based on the observed centroids.

`-pad num_pad_images`

This overrides the number of images (default=0) that will be added *before and after* the predicted number of images that a reflection needs for its 3D profile. The extra images can be useful for background determination and to allow for crystal slippage. The actual number of images used for a reflection profile depends on the effective mosaic spread and reflection position (Lorentz factor).

`-window fast_size slow_size`

This overrides the nominal size of the window of pixels (default=21 by 21 pixels) in the image that will be accumulated in the 3D reflection profile. The actual size for any given reflection may be adjusted internally for spatial distortion and oblique incidence. The sizes need to be large enough to include a local background. The sizes can be so large that neighboring reflections are included in the area. These neighbors will be excluded automatically from the background determination.

`-profsize profile_size1 profile_size2`

This overrides the size of the Kabsch-style 3D profile (default=0.5, 0.3) for reflections. The first parameter describes the degrees subtended by the reflection profile in a plane perpendicular to the scattered beam wavevector, while the second parameter describes the degrees subtended by the reflection profile along a vector from the origin of the Ewald sphere and bisecting the reciprocal lattice vector. The first parameter is mainly a function of the reflection size on the detector and the crystal to detector distance. It can be approximated by  $\arctan(\text{window size} / \text{crystal-to-detector distance})$ . The second parameter is mainly a function of the

effective mosaic spread of the crystal, but in no case should be less than the rotation width of a single image.

`-nonunf file_with_nonunf_info_in_header`

Specifies an alternate file to retrieve non-uniformity information from. The default info comes from the `-scan` file. This option overrides the information in there. Note that this does not *directly* specify the non-uniformity file to use. The *header* of this file specifies the information to use which in turn may point to another file (which can be the same file).

## 4.2 Examples

```
dtintegrate -scan refine.head -seq 0 10
dtprofit
```

The above command line reads image `refine.head` for scan information and overrides the image sequence range found in that image with a new range of 0 to 10. Any non-uniformity correction information is found in the header of image `refine.head`. Results are written to the file `reflnprofiles.profit` which are then read by the `dtprofit` program to create the file `dtprofit.ref`.

```
dtintegrate -scan refine.head -reso 3.5 50 -pad 2 -batch 3 -window 25 25
dtprofit
```

The first command reads image `refine.head` for scan information and overrides the default resolution value of the entire detector surface with a resolution range of 3.5 to 50 Ångstroms. Each reflection 3D profile is padded by 2 images and is 25 by 25 pixels in size. Refinement is carried out after every 3 images.

## 4.3 Results

The command:

```
dtintegrate -scan newrefine.head -reso 10. 3.5 -pad 2 \
            -window 25 25 -batch 3 -seq 1 44
dtprofit
```

produced the following output which is explained below.

```
1  File newrefine.head successfully opened.
2  Cnonunf::nInitValues called
3  D0_NONUNF_TYPE: >>None<<
4  Command line string: >>-reso<<
```

```
5 Command line string: >>-pad<<
6 Command line string: >>-window<<
7 Command line string: >>-batch<<
8 Command line string: >>-seq<<
9 Integrate object listing:
10 Pad: 2
11 Window size: 25, 25
12 Profile size: 0.5, 0.3
13 Scan seqrage: 1, 44
14 Images/batch: 3
15 Resolution: 10, 3.5
16 dtintegrate: 3D method used
17 Cnonunf::nInitValues called
18 DO_NONUNF_TYPE: >>None<<

19 File /u6/jwp/DATA/NUCL4/nucl4001.osc successfully opened.
20 Header of file /u6/jwp/DATA/NUCL4/nucl4001.osc successfully read.
21 Image is 1900 by 1900 pixels.
22 Data_type in header is unsigned short int.
23 Compression_type is RAXIS.
24 Byte_order is big_endian.

25 IMAGE #: 0 Name: /u6/jwp/DATA/NUCL4/nucl4001.osc

26 Rotation list:
27 Start: 24
28 End: 24.25
29 Increment: 0.25
30 Time: 60
31 Predicting reflns...
32 Done predicting reflns...
33 Processing reflns ...
34 there are 597 in the list.
35 Number new, active, full, donew, dones:
36 206, 206, 0, 0, 0
37 Reflections deleted: 391

38 File /u6/jwp/DATA/NUCL4/nucl4002.osc successfully opened.
39 Header of file /u6/jwp/DATA/NUCL4/nucl4002.osc successfully read.
40 Image is 1900 by 1900 pixels.
41 Data_type in header is unsigned short int.
42 Compression_type is RAXIS.
43 Byte_order is big_endian.

44 IMAGE #: 1 Name: /u6/jwp/DATA/NUCL4/nucl4002.osc

45 Rotation list:
46 Start: 24.25
47 End: 24.5
48 Increment: 0.25
49 Time: 60
50 Predicting reflns...
51 Done predicting reflns...
52 Processing reflns ...
53 there are 803 in the list.
54 Number new, active, full, donew, dones:
55 72, 278, 0, 0, 0
56 Reflections deleted: 525

57 File /u6/jwp/DATA/NUCL4/nucl4003.osc successfully opened.
58 Header of file /u6/jwp/DATA/NUCL4/nucl4003.osc successfully read.
59 Image is 1900 by 1900 pixels.
60 Data_type in header is unsigned short int.
61 Compression_type is RAXIS.
62 Byte_order is big_endian.

63 IMAGE #: 2 Name: /u6/jwp/DATA/NUCL4/nucl4003.osc
```

```
64 Rotation list:
65 Start: 24.5
66 End: 24.75
67 Increment: 0.25
68 Time: 60
69 Predicting reflns...
70 Done predicting reflns...
71 Processing reflns ...
72 there are 858 in the list.
73 Number new, active, full, donew, dones:
74 56, 334, 0, 0, 0
75 Reflections deleted: 524
76 ... end of batch ...

77 File /u6/jwp/DATA/NUCL4/nucl4004.osc successfully opened.
78 Header of file /u6/jwp/DATA/NUCL4/nucl4004.osc successfully read.
79 Image is 1900 by 1900 pixels.
80 Data_type in header is unsigned short int.
81 Compression_type is RAXIS.
82 Byte_order is big_endian.

83 IMAGE #: 3 Name: /u6/jwp/DATA/NUCL4/nucl4004.osc

84 Rotation list:
85 Start: 24.75
86 End: 25
87 Increment: 0.25
88 Time: 60
89 Predicting reflns...
90 Done predicting reflns...
91 Processing reflns ...
92 there are 1139 in the list.
93 Number new, active, full, donew, dones:
94 80, 414, 0, 0, 0
95 Integrate: start...
96 Strong: 0, Weak: 0
97 Integrate: ...done.
98 Reflections deleted: 725

99 File /u6/jwp/DATA/NUCL4/nucl4005.osc successfully opened.
100 Header of file /u6/jwp/DATA/NUCL4/nucl4005.osc successfully read.
101 Image is 1900 by 1900 pixels.
102 Data_type in header is unsigned short int.
103 Compression_type is RAXIS.
104 Byte_order is big_endian.

105 IMAGE #: 4 Name: /u6/jwp/DATA/NUCL4/nucl4005.osc

106 Rotation list:
107 Start: 25
108 End: 25.25
109 Increment: 0.25
110 Time: 60
111 Predicting reflns...
112 Done predicting reflns...
113 Processing reflns ...
114 there are 1018 in the list.
115 Number new, active, full, donew, dones:
116 70, 484, 20, 0, 20
117 Reflections deleted: 534

118 File /u6/jwp/DATA/NUCL4/nucl4006.osc successfully opened.
119 Header of file /u6/jwp/DATA/NUCL4/nucl4006.osc successfully read.
120 Image is 1900 by 1900 pixels.
121 Data_type in header is unsigned short int.
122 Compression_type is RAXIS.
123 Byte_order is big_endian.
```

```
124 IMAGE #: 5 Name: /u6/jwp/DATA/NUCL4/nucl4006.osc
125 Rotation list:
126 Start: 25.25
127 End: 25.5
128 Increment: 0.25
129 Time: 60
130 Predicting reflns...
131 Done predicting reflns...
132 Processing reflns ...
133 there are 1086 in the list.
134 Number new, active, full, donew, dones:
135 70, 534, 40, 0, 60
136 Reflections deleted: 532
137 ... end of batch ...
138 Refine listing:

139 Crystal listing:

140 Unit cell lengths: 73.551, 68.265, 104.777
141 Unit cell angles: 90, 90, 90
142 Unit cell volume: 526081

143 Spacegroup number: 16
144 name: P222
145 Num. equiv. posns: 4
146 Equival. position 1:
147     0,     0     0
148     1,     0     0
149     0,     1     0

150 Equival. position 2:
151     0,     0     0
152    -1,     0     0
153     0,     1     0

154 Equival. position 3:
155     0,     0     0
156    -1,     0     0
157     0,    -1     0

158 Equival. position 4:
159     0,     0     0
160     1,     0     0
161     0,    -1     0

162 Orientation angles: 172.744, -29.3, 15.008
163 Mosaic Spread: 0.6
164 Description: Test of d*TREK modules

165 CRYSTAL_Goniometer listing:

166 Description: A 3-circle eulerian crystal goniostat
167 Number of values: 3

168 Value 0 name: Omega
169 units: deg
170 Vector: 1, 0, 0
171 Datum value: -24
172 Current value: Unknown

173 Value 1 name: Chi
174 units: deg
175 Vector: 0, 0, 1
176 Datum value: 0
177 Current value: Unknown

178 Value 2 name: Phi
```

```
179 units: deg
180 Vector: 1, 0, 0
181 Datum value: 0
182 Current value: Unknown

183 Source listing:

184 Single wavelength: 1.54178
185 Direction vector: 0, 0, 1
186 Spectral Dispersion: 0.0002, 0.0002
187 Crossfire: 0, 0, 0, 0
188 Polarization: 0.5, 1, 0, 0
189 Size: 0, 0, 0, 0
190 Intensity: 0
191 Refinement flags: Not yet available

192 D0_Detector listing:

193 Pixel dimensions: 1900, 1900
194 Nominal size in mm: 193.23, 199.5
195 Description: R-AXIS IIC at MSC
196 Rotation: 0, 0, 0
197 Translation: -0.142, 0.197, -96.346
198 Fast pixel vector: 0, 1, 0
199 Slow pixel vector: -1, 0, 0

200 D0_Goniometer listing:

201 Description: A simple detector goniostat
202 Number of values: 6

203 Value 0 name: RotZ
204 units: deg
205 Vector: 0, 0, 1
206 Datum value: 1.234
207 Current value: Unknown

208 Value 1 name: RotX/Swing
209 units: deg
210 Vector: 1, 0, 0
211 Datum value: 0
212 Current value: Unknown

213 Value 2 name: RotY
214 units: deg
215 Vector: 0, 1, 0
216 Datum value: 0
217 Current value: Unknown

218 Value 3 name: TransX
219 units: mm
220 Vector: 1, 0, 0
221 Datum value: -0.142
222 Current value: Unknown

223 Value 4 name: TransY
224 units: mm
225 Vector: 0, 1, 0
226 Datum value: 0.197
227 Current value: Unknown

228 Value 5 name: TransZ/Dist
229 units: mm
230 Vector: 0, 0, -1
231 Datum value: 96.346
232 Current value: Unknown
```

```
233 DetResolution min: 154178
234 DetResolution max: 1.65533

235 Internal loop:      0
236 Reflections accepted: 60
237 Reflections rejected: 0

238 Internal loop:      1
239 Reflections accepted: 60
240 Reflections rejected: 0

241 Internal loop:      2
242 Reflections accepted: 60
243 Reflections rejected: 0

244 Internal loop:      3
245 Reflections accepted: 60
246 Reflections rejected: 0

247 Internal loop:      4
248 Reflections accepted: 60
249 Reflections rejected: 0

250 Refinement residuals
251 rmsResid (mm) = 0.130414
252 rmsResid (Deg) = 0.0344216

253 Crystal
254 a, b, c:      73.395      68.225      104.852
255 Sigmas:      0.0646      0.0110      0.0252
256 Shifts:      0.0000      0.0000      0.0001
257 alpha, beta, gamma:  90.000      90.000      90.000
258 Sigmas:      0.0000      0.0000      0.0000
259 Shifts:      0.0000      0.0000      0.0000
260 Rot1, Rot2, Rot3:  172.748      -29.297      14.914
261 Sigmas:      0.0036      0.0100      0.0726
262 Shifts:      0.0000      -0.0002      -0.0249
263 Mosaicity:      0.600
264 Sigma:      0.0000
265 Shift:      0.0000

266 Detector: 0
267 Trans1, Trans2, Trans3:  -0.008      0.158      96.996
268 Sigmas:      0.1148      0.0459      0.1999
269 Shifts:      0.0000      0.0000      0.0210
270 Rot1, Rot2, Rot3:  0.848      0.000      0.000
271 Sigmas:      0.0837      0.0000      0.0000
272 Shifts:      -0.0125      0.0000      0.0000

273 Source
274 Rot1, Rot2, Wavelength:  0.000      0.000      1.5418
275 Sigmas:      0.0000      0.0000      0.0000
276 Shifts:      0.0000      0.0000      0.0000

277 File intrefine.head successfully opened.

278 File /u6/jwp/DATA/NUCL4/nucl4007.osc successfully opened.
279 Header of file /u6/jwp/DATA/NUCL4/nucl4007.osc successfully read.
280 Image is 1900 by 1900 pixels.
281 Data_type in header is unsigned short int.
282 Compression_type is RAXIS.
283 Byte_order is big_endian.

284 IMAGE #: 6 Name: /u6/jwp/DATA/NUCL4/nucl4007.osc

285 Rotation list:
286 Start: 25.5
287 End: 25.75
```

```
288 Increment: 0.25
289 Time: 60
290 Predicting reflns...
291 Done predicting reflns...
292 Processing reflns ...
293 there are 1369 in the list.
294 Number new, active, full, donew, dones:
295 74, 568, 58, 0, 118
296 Integrate: start...
297 Strong: 118, Weak: 0
298 Integrate: ...done.
299 Reflections deleted: 859

300 File /u6/jwp/DATA/NUCL4/nucl4008.osc successfully opened.
301 Header of file /u6/jwp/DATA/NUCL4/nucl4008.osc successfully read.
302 Image is 1900 by 1900 pixels.
303 Data_type in header is unsigned short int.
304 Compression_type is RAXIS.
305 Byte_order is big_endian.

306 IMAGE #: 7 Name: /u6/jwp/DATA/NUCL4/nucl4008.osc

307 Rotation list:
308 Start: 25.75
309 End: 26
310 Increment: 0.25
311 Time: 60
312 Predicting reflns...
313 Done predicting reflns...
314 Processing reflns ...
315 there are 1110 in the list.
316 Number new, active, full, donew, dones:
317 63, 573, 64, 0, 64
318 Reflections deleted: 537

319 File /u6/jwp/DATA/NUCL4/nucl4009.osc successfully opened.

320 [... many lines deleted ...]

321 IMAGE #: 42 Name: /u6/jwp/DATA/NUCL4/nucl4043.osc

322 Rotation list:
323 Start: 34.5
324 End: 34.75
325 Increment: 0.25
326 Time: 60
327 Predicting reflns...
328 Done predicting reflns...
329 Processing reflns ...
330 there are 1502 in the list.
331 Number new, active, full, donew, dones:
332 69, 603, 71, 0, 215
333 Integrate: start...
334 Strong: 215, Weak: 0
335 Integrate: ...done.
336 Reflections deleted: 970

337 File /u6/jwp/DATA/NUCL4/nucl4044.osc successfully opened.
338 Header of file /u6/jwp/DATA/NUCL4/nucl4044.osc successfully read.
339 Image is 1900 by 1900 pixels.
340 Data_type in header is unsigned short int.
341 Compression_type is RAXIS.
342 Byte_order is big_endian.

343 IMAGE #: 43 Name: /u6/jwp/DATA/NUCL4/nucl4044.osc

344 Rotation list:
345 Start: 34.75
```

```
346 End: 35
347 Increment: 0.25
348 Time: 60
349 Predicting reflns...
350 Done predicting reflns...
351 Processing reflns ...
352 there are 1150 in the list.
353 Number new, active, full, donew, dones:
354 77, 609, 68, 0, 68
355 Reflections deleted: 541
```

## Explanation

- Line 1            `dtintegrate` is started and file `newrefine.head` is read.
- Lines 2-3        The non-uniformity object is initialized based on the information in `newrefine.head`. In this case, there is no non-uniformity to be applied.
- Lines 4-8        Command line options are parsed.
- Lines 9-18       Parameters used by `dtintegrate` are listed. The values are either the default values or the ones overridden on the command line.
- Lines 19-24     The first image in the scan is read in and information from its header is displayed.
- Line 25          The first image is numbered 0 and its name is listed.
- Lines 26-30     The rotation information from the image is listed. This includes the rotation start, end, increment and time.
- Lines 31-34     Reflections which appear on the image are predicted.
- Lines 35-37     Reflections which appear on the image are processed in a first pass. There were 206 new 3D arrays (shoeboxes) allocated and data from the image moved into them (`new`). There were a total of 206 active reflections (`active`) open in this image. No shoeboxes were filled. (`full`). Since there were no full shoeboxes, none were processed either as weak (`donew`) or strong (`dones`) reflections. Just remember that if you have all weak reflections, then you will not be able to do profile fitting, as you have no strong reflections to compute reference profiles from.
- Line 55          There were 56 new reflections and 334 active reflections in image 2.

Line 116      Finally some shoeboxes were filled up. There were 20 full shoeboxes and there were all preprocessed as strong.

Lines 138-276 This is the refinement loop. Even though `-batch 3` was specified on the command line, there were no reflections available for refinement until the second batch of images was finished.

Consult the `tCrefine` documentation for an explanation of all of these lines.

Lines 251-252 The root mean square difference between predicted and observed reflections centroids is listed.

Line 277      `dtintegrate` always writes the refinement results to the file `intrefine.head`.

Line 278      Integration resumes with reading in image 6 with sequence number 7.

Line 319      Many lines were deleted because they are repetitive.

Line 342      The last image with sequence number 44 is read in and processed.

The file `dtintegrate.ref` with the background-peak-background integration of the reflections was written. It has a standard `d*TREK` reflection file format. Below is a truncated listing of the file produced by the above example.

```
1  5 14 0
2  nH
3  nK
4  nL
5  nDetector_number
6  nNonunf_flag
7  fIntensity
8  fSigmaI
9  fCalc_pixel0
10 fCalc_pixell
11 fCalc_rot_mid
12 fCalc_rot_width
13 fCalc_polarz
14 fCalc_lorentz
15 fCalc_oblique
16 fResolution
17 fObs_pixel0
18 fObs_pixell
19 fObs_rot_mid
20 fObs_rot_width
21      3 -17 -7 0 0 85954.8 602.587 1359.54 1010 24.5148 0.614138 0.924924 2.58067
    1.08823 3.82877 1365.76 1010.82 24.5198 -999
22     -3 10 -3 0 0 364503 435.747 734.539 974.991 24.9057 0.607173 0.972394 4.25585
    1.02933 6.45861 731.83 975.897 24.9105 -999
```

```
23      3  18  -1  0  0  660363  780.311  565.905  823.491  24.5262  0.63487  0.926624  2.6104
1.09269  3.74323  559.694  821.784  24.5122  -999
24      9 -16  -2  0  0  41327.2  634.775  1364.66  863.393  24.8027  0.613915  0.922555  2.54091
1.09131  3.76893  1370.44  862.447  24.7832  -999
25      7 -13  -1  0  0  22344.9  387.538  1278.05  881.164  24.3121  0.612952  0.949135  3.13527
1.05716  4.68836  1282.75  880.649  24.3156  -999
```

The output of the dtprofit program follows.

```
1  !***ERROR in DSKBR, reading error.
2  !***ERROR in DSKBR, reading error.
3  Number of reflection profiles read in: 2625
4  Sum (Pi**2/Vi) = 1001.34
5  File D0000.krefer successfully opened.
6  Success writing file D0000.krefer!

7  Sigma Range   Avg.Corr.   Num.Refln.
8  0... 3        -0.18      26
9  3... 6         0.17      11
10 6... 9         0.164     13
11 9...12        0.218     14
12 12...15       0.255     14
13 15...18       0.301     17
14 18...21       0.376     15
15 21...24       0.315     11
16 24...27       0.324      8
17 27...30       0.842    2496

18 Corr.Range   Num.Refln.
19 0...0.1      57
20 0.1...0.2   28
21 0.2...0.3   39
22 0.3...0.4   53
23 0.4...0.5   63
24 0.5...0.6  135
25 0.6...0.7  199
26 0.7...0.8  253
27 0.8...0.9  339
28 0.9... 1   1459
```

Line 1        The file reflnprofiles.profit generated by dtintegrate was read. Reflection profiles were processed until an end-of-file was reached. Reference profiles were calculated.

Line 2        The profiles file was read a second time until an end-of-file was reached. Observed profiles were analyzed with the Kabsch method.

Line 3        There were 2625 reflection profiles in the input file.

Line 4        The sum of the square of the variance weighted reference profiles is 1001.34. This number has no special meaning.

Lines 5-6     A d\*TREK image file was written with the overall reference profile in it. This image file can be viewed with dtdisplay as a diagnostic tool.

Lines 7-17 The correlation of the observed profiles to the reference profile in 10 different ranges of  $I/\sigma I$  is given. The lower ranges should not have a very high correlation, while the higher ranges should have a correlation above 0.8. These correlations confirm that the assumption that the reflections have the same profiles.

Lines 18-28 The number of reflections in 10 different correlation ranges is given. This also confirms the underlying assumptions of profile analysis.

The reflection results of `dtprofit` are written to the file `dtprofit.ref` which is a standard `d*TREK` reflection file format. Below is a truncated example of this file. The fields `fIntensity` and `fSigmaI` contain the profile analyzed reflection intensity. The fields `fOtherInt` and `fOtherSig` contain the background-peak-background intensity estimated by `dtintegrate`. All intensities and sigmas have been corrected for both the Lorentz and polarization factors which is also included in the field `fLPfactor` as the multiplier.

```
1  4 12 0
2  nH
3  nK
4  nL
5  nBadFlag
6  fIntensity
7  fSigmaI
8  fOtherInt
9  fOtherSig
10 fObs_pixel0
11 fObs_pixel1
12 fObs_rot_mid
13 fCalc_pixel0
14 fCalc_pixel1
15 fCalc_rot_mid
16 fLPfactor
17 fCorrelation
18   3 -17 -7 0 86174.8 274.128 85954.8 602.587 1365.76 1010.82 24.5198 1359.54
19 1010 24.5148 0.418949 1.00256
20 -3 10 -3 0 365861 322.356 364503 435.747 731.83 975.897 24.9105 734.539
21 974.991 24.9057 0.241641 1.00373
22   3 18 -1 0 644888 574.647 660363 780.311 559.694 821.784 24.5122 565.905
23 823.491 24.5262 0.413418 0.976566
24   9 -16 -2 0 35310.3 229.748 41327.2 634.775 1370.44 862.447 24.7832 1364.66
25 863.393 24.8027 0.426598 0.854408
26   7 -13 -1 0 21440.3 151.118 22344.9 387.538 1282.75 880.649 24.3156 1278.05
27 881.164 24.3121 0.336045 0.959517
```

#### 4.4 Overcoming difficulties

Since `dtintegrate` gets all of its input information from an image header, difficulties may arise when the header is incomplete. Be sure that the detector spatial distortion information and non-uniformity information is correct and the files exist in the current directory or the filenames in the header include the path or directory.



```
float a3fScatt[3];
float fRecipLP;
} tagProfitRec;
```

followed by the floating point profile (grid1 \* grid2 \* grid3 \* 4 bytes) and the floating point variance in the profile (grid1 \* grid2 \* grid3 \* 4 bytes).